

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	<u>Salt 1</u> <u>Preliminary Test</u>		
1)	Note the colour of the given salt.	Colourless	Absence of Cu^{2+} , Ni^{2+} , Mn^{2+} , Fe^{2+}
2)	Note the smell of the given salt.	No characteristic smell	Absence of NH_4^+ , CH_3COO^-
3)	Solubility in dil. HCl	Soluble	Absence of Pb^{2+}
4)	Flame test A pinch of salt is made into paste with conc. HCl and the paste is introduced to the flame by using a glass rod.	Brick red flame	Presence of Ca^{2+}

5) Ash test

The test is as follows:

No characteristic

Absence of Al^{3+} , Zn^{2+} , Mg^{2+}

A pinch of salt is made into paste with conc. HCl and the paste is introduced to the flame by using a glass rod.

Back near flame

Presence of

5) Ash test

The test solution is mixed with few drops of conc. HNO_3 and a little cobalt nitrate solution. A rolled filter paper dipped in to it is burnt to ash.

No characteristic colour in the ash


Presence of Al^{3+} , Zn^{2+} , Mg^{2+}

Analysis of anion

1) To the test solution add dil. HCl.

Gives effervescence of colourless, odourless, gas which turns lime water milky.

Presence of carbonate (CO_3^{2-})



SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
2)	The test solution is rubbed with dil. H_2SO_4	No smell of vinegar	Absence of acetate (CH_3COO^-)
3)	To a little of the salt solution about 1ml of conc. H_2SO_4 is added & heated.	No dense white gas.	Absence of chloride (Cl^-)
4)	A little of the salt soln is heated with conc. H_2SO_4 & a paper ball is added	No reddish brown fumes	Absence of nitrate (NO_3^-)
5)	To the test solution add $BaCl_2$ solution.	No white precipitate.	Absence of sulphate (SO_4^{2-})



BaCl₂ solution.

sulphate (SO₄²⁻)


Confirmatory Test

Carbonate

- 1) To the test solution add white precipitate
BaCl₂ solution which is soluble
in dil. HCl. Presence of CO₃²⁻
confirmed



- 2) A little of the solid salt is heated with dil. H₂SO₄ Colourless, odourless gas which turns lime water milky. Presence of CO₃²⁻
confirmed.



SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	<u>Analysis of cations Intergroup separation</u>		
1)	To the test solution add dil Na_2CO_3	white precipitate	Absence of I group (NH_4^+)
2)	To the test solution add dil HCl	No white precipitate	Absence of II group (Pb^{2+})
3)	To the test solution add dil HCl & pass H_2S gas.	No black precipitate	Absence of III group (Cu^{2+})
4)	To the test solution add dil HCl , NH_4Cl & NH_4OH pass H_2S gas.	No white gelatinous precipitate	Absence of III group (Al^{3+} , Fe^{2+})
5)	To the test soln add dil HCl , NH_4Cl & NH_4OH pass H_2S gas.	No precipitate	Absence of IV group (Zn^{2+} , Mn^{2+})


6)	To the test soln add dil HCl , NH_4Cl , NH_4OH	white precipitate	Presence of V group (Ca^{2+} , Ba^{2+} , Sr^{2+})
----	---	-------------------	--

4) To the test solution add dil. HCl, NH_4Cl & NH_4OH pass H_2S gas.	No white gelatinous precipitate.	Absence of III group (Al^{3+} , Fe^{2+})
5) To the test soln add dil. HCl, NH_4Cl & NH_4OH pass H_2S gas.	No precipitate	Absence of IV group (Zn^{2+} , Mn^{2+})


6) To the test soln add dil. HCl, NH_4Cl , NH_4OH and $(\text{NH}_4)_2\text{CO}_3$	white precipitate	Presence of V group (Ca^{2+} , Ba^{2+} , Sr^{2+})
7) In the test soln add dil. HCl, NH_4Cl , NH_4OH & Na_2HPO_4	No white precipitate	Absence of VI group (Mg^{2+})

Confirmatory test for fifth group (Ca^{2+} , Ba^{2+} , Sr^{2+})

1) To the test soln acetic acid & potassium chromate soln is added	Yellow solution	Presence of Ca^{2+}
--	-----------------	------------------------------



2) Ammonium oxalate solution and NH_4OH soln is added to the test soln. To the above precipitate acetic acid is added & heated.	white precipitate insoluble in acetic acid	Presence of Ca^{2+}
---	--	------------------------------



SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	<u>Salt 2</u> <u>Preliminary Test</u>		
1.	Noted the colour of the given salt	Colourless	Absence of Cu^{2+} , Ni^{2+} , Mn^{2+} , Fe^{2+}
2.	Noted the smell of the given salt	No characteristic smell	Absence of NH_4^+ , CH_3COO^-
3.	Solubility in dil. HCl.	Insoluble	Presence of Pb^{2+}
4.	Flame test A pinch of salt solution is made to paste with conc. HCl and it is introduced to flame by using a glass rod.	No characteristic colour on flame.	Absence of Ca^{2+} , Sr^{2+} , Ba^{2+}

5.	Asb test The test solution is	No characteristic colour on ash	Absence of Al^{3+} , Zn^{2+} , Mg^{2+}
----	----------------------------------	---------------------------------	---

4. Flame test

A pinch of salt solution is made to paste with conc. HCl and it is introduced to flame by using a glass rod.

No characteristic colour on flame.

Absence of Ca^{2+} , Sr^{2+} & Ba^{2+}

5. Ash test

The test solution is mixed with few drops of conc. HNO_3 and a little cobalt nitrate solution. A rolled filter paper dipped in to it is burnt to ash.

No characteristic colour on ash

Absence of Al^{3+} , Zn^{2+} , Mg^{2+}


Analysis of anion

1. To the test solution add dil. HCl. No brisk effervescence.

Absence of carbonate (CO_3^{2-}).

2. To the test solution add dil. H_2SO_4 . Smell of vinegar.

Presence of acetate (CH_3COO^-).



3. To the little of the salt solution is heated with conc. H_2SO_4 and a paper ball. No reddish brown fumes.

No reddish brown fumes

Absence of nitrate (NO_3^-)

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	is added.		
4.	To the little of the salt solution about 1ml of conc. H_2SO_4 is added & heated.	No dense white gas.	Absence of chloride (Cl^-).
5.	To the test solution add $BaCl_2$ solution.	No white precipitate	Absence of sulphate (SO_4^{2-})
	<u>Confirmatory test</u> <u>Acetate</u>		
1.	To the test solution add neutral ferric chloride solution.	Deep red colouration	Presence of acetate.
	$3CH_3COOH + FeCl_3 \rightarrow (CH_3COO)_3Fe + 3HCl$		


2.	To the test solution add ethanol and 1ml conc. H_2SO_4	Pleasant fruity.	Presence of acetate.
----	--	------------------	----------------------

Confirmatory test Acetate

1. To the test solution add neutral ferric chloride solution.

Deep red colouration


Presence of acetate.



2. To the 1st solution add ethanol and 1ml conc. H_2SO_4 . Heat the mixture and pour the contents into a beaker containing water.

Pleasant fruity.

Presence of acetate.



Analysis of cation intergroup separation

1. To the test solution add dil. Na_2CO_3 .

white precipitate

Absence of O group (NH_4^+)

2. To the test solution add dil. HCl .

white precipitate.

Presence of I group (Pb^{2+})

3. To the test solution add dil. HCl and pass H_2S gas.

No Black precipitate.

Absence of II group (Cu^{2+})

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
4.	To the test solution add dil. HCl, NH_4Cl and NH_4OH .	No white precipitate.	Absence of III group (Al^{3+} , Fe^{2+})
5.	To the test solution add dil. HCl, NH_4Cl and NH_4OH pass H_2S gas.	No white & black precipitate.	Absence of IV group (Ca^{2+} , Ba^{2+} , Sr^{2+})
6.	To the test solution add dil. HCl, NH_4Cl , NH_4OH and $(\text{NH}_4)_2\text{CO}_3$.	No white precipitate.	Absence of V group (Zn^{2+} , Mn^{2+})
7.	To the test solution add dil. HCl, NH_4Cl , NH_4OH and Na_2HPO_4 .	No white precipitate.	Absence of VI group (Mg^{2+})
Confirmatory test for lead ion (Pb^{2+})			

1.) To the test solution is acidified with acetic	Yellow precipitate	Presence of Pb^{2+}
---	--------------------	------------------------------

6. To the test solution add dil HCl, NH_4Cl , NH_4OH and $(\text{NH}_4)_2\text{CO}_3$.	No white precipitate.	Absence of IV group (Zn^{2+} , Mn^{2+})
7. To the test solution add dil HCl, NH_4Cl , NH_4OH and Na_2HPO_4 .	No white precipitate.	Absence of VI group (Mg^{2+})
Confirmatory test for lead ion (Pb^{2+})		

1) To the test solution is acidified with acetic acid & potassium chromate is added	Yellow precipitate	Presence of Pb^{2+}
$\text{Pb}^{2+} + \text{K}_2\text{CrO}_4 \rightarrow \text{PbCrO}_4 + 2\text{K}^+$		
2) Potassium iodide solution is added to the test solution	Yellow precipitate	Presence of lead
$\text{Pb}^{2+} + \text{K}_2\text{CrO}_4 \rightarrow \text{PbI}_2 \downarrow + 2\text{K}^+$		
The above precipitate is boiled with water & cooled.	Precipitate dissolves on boiling. Golden yellow sparkles is found on cooling.	Presence of lead confirmed.
<u>Result</u> The given salt contain lead as cation & acetate as anion. Therefore the salt is lead acetate <u>lead acetate</u> . 25/05/20		

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	<u>Salt 3</u>		
	<u>Preliminary Test</u>		
1.	Note the colour of the given salt.	Colourless	Absence of Cu^{2+} , Ni^{2+} , Mn^{2+} , Fe^{2+}
2.	Note the smell of the given salt.	Ammoniacal smell	May be the presence of NH_4^+
3.	Solubility in dil. HCl	Soluble	Absence of Pb^{2+}
4.	Flame test A pinch of salt is made into paste with conc. HCl and the paste is introduced to the flame by using a glass rod.	No characteristic colour on the flame.	Absence of Ca^{2+} , Sr^{2+} , Ba^{2+}

5.	Ash test The solution is mixed with a few drops of conc. HCl and the paste is introduced to the flame by using a glass rod.	No characteristic colour on the flame.	Absence of Al^{3+} , Zn^{2+} , Mg^{2+}
----	--	--	---

into paste with conc. HCl and the paste is introduced to the flame by using a glass rod.

flame.

5. Ash test

The solution is mixed with few drops of conc. HNO_3 and a little cobalt nitrate soln. A rolled filter paper dipped in it is burnt to ash.

No characteristic colour on the ash.

Absence of Al^{3+} , Zn^{2+} , Mg^{2+}

Analysis of Anion

1. To the test solution add dil. HCl.

No brisk effervescence.

Absence of carbonate (CO_3^{2-}).

2. The test solution is rubbed with dil. H_2SO_4 .

No vinegar like smell.

Absence of acetate (CH_3COO^-).

3. To a little of the salt soln about 1ml of conc. H_2SO_4 is added and heated.

Dense white gas which fumes with NH_4OH solution taken on a glass rod.

Presence of chloride (Cl^-).

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	$MCl + H_2SO_4 \rightarrow MH_2SO_4 + HCl$ $MCl + NH_4OH \rightarrow NH_4Cl + H_2O$		
4.	A little of the salt soln is heated with conc H_2SO_4 and a paper ball is added.	No Reddish brown fumes	Absence of nitrate (NO_3^-)
5.	To the test solution add bally soln	No white precipitate	Absence of sulphate (SO_4^{2-})
6.	<u>Confirmatory Test.</u>		
	<u>Chloride</u>		
1.	To the test solution is heated with conc H_2SO_4 and a pinch of manganese dioxide (MnO_2).	Greenish yellow gas with irritating smell.	Presence of chloride.



2.	The test solution is mixed with dil. H_2O	White curdy precipitate soluble	Presence of chlorine combined
----	---	---------------------------------	-------------------------------

5. To the test solution add BaCl_2 soln

No white precipitate

Absence of sulphate (SO_4^{2-})

Confirmatory Test

Chloride

1. To the test solution is heated with conc H_2SO_4 and a pinch of manganese dioxide

Greenish yellow gas with irritating smell



Presence of chloride



2. The test solution is mixed with dil. HNO_3 and add AgNO_3 solution

White curdy precipitate soluble in excess of NH_4OH

Presence of chlorine confirmed.



Analysis of cation (Inter-group separation)

1) To the test soln add dil. Na_2CO_3

No precipitate

Presence of I group (NH_4^+)

2) To the test soln add dil. HCl

No white precipitate

Absence of II group (Pb^{2+})

3) To the test soln add dil. HCl , NH_4Cl , NH_4OH

No black precipitate

Absence of III group (Cu^{2+})

4) To the test soln add dil. HCl , NH_4Cl , NH_4OH

No white gelatinous precipitate

Absence of IV group (Al^{3+} , Fe^{2+})

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
5)	To the test soln add dil. HCl, NH_4Cl and NH_4OH gas H_2S .	No white precipitate.	Absence of IV group (Mn^{2+} , Zn^{2+})
6)	To the test soln add dil. HCl, NH_4Cl and NH_4OH and $(\text{NH}_4)_2\text{CO}_3$.	No white precipitate.	Absence of V group (Ca^{2+} , Ba^{2+} , Sr^{2+})
7)	To the test solution add dil. HCl, NH_4Cl , NH_4OH and Na_2HPO_4 .	No white precipitate.	Absence of VI group (Mg^{2+}).
Confirmatory test for ammonium ion			
1)	NaOH soln is added the test solution is excess	Colourless gas with smell of ammonia. The gas gives white fumes	Absence of NH_4^+

when a glass rod is held in dil.

and Na_2HPO_4

precipitate


(Hg^{2+})

Confirmatory test for ammonium ion

1) NaOH soln is added the test solution is concn

colorless gas with small amount


when a glass rod dipped in dil. HCl is brought near the mouth



2) Nessler's reagent is added to the test soln.

Brown precipitate

Presence of NH_4^+



Result

The result salt contain ammonium as cation and chloride as anion.

25/10/20

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	<u>Salt 4.</u>		
1.	Note the colour of the given salt.	Colourless	Absence of Cu^{2+} , Ni^{2+} , Mn^{2+} , Fe^{2+}
2.	Note the smell of the given salt.	No characteristic smell	Absence of NH_4^+ , CH_3COO^-
3.	Solubility in dil. HCl	Soluble	Absence of Pb^{2+}
4.	Flame test A pinch of salt is mixed with conc. HCl & colour on the flame Ca^{2+} , Ba^{2+} , the paste is introduced as to the flame by using a glass rod.	NO characteristic colour	Absence of Ca^{2+} , Ba^{2+}
5.	Ash test The test solution is	Green tinted ash	Presence of Zinc

mixed with few drops of conc. HNO_3 and a little cobalt nitrate solution
A rolled filter paper

4. Flame test

A pinch of salt, besides NO characteristic Absence of Ca^{2+} , into paste with conc. HCl & colour on the flame Sr^{2+} , Ba^{2+} , the paste is introduced as to the flame by using a glass rod.

5. Ash test

The test solution is

Green tinted ash

Absence of Zinc

mixed with few drops of conc. HNO_3 and a little cobalt nitrate solution. A rolled filter paper dipped in to it is burnt to ash.

Analysis of anion

1. To the test solution add dil. HCl

No brisk effervescence.

Absence of carbonate (CO_3^{2-})

2. To the test solution is mixed with dil. H_2SO_4

No smell of vinegar.

Absence of acetate (CH_3COO^-)

3. To a little of the salt soln about 1ml of conc. H_2SO_4 is added & heated.

No dense white gas

Absence of chloride (Cl^-)

4. A little of the salt soln is heated with conc. H_2SO_4 and paper ball is added.


No reddish brown fumes.

Absence of nitrate (NO_3^-)

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
5.	To the test solution add BaCl_2 solution $\text{M}_2\text{SO}_4 + \text{BaCl}_2 \rightarrow \text{BaSO}_4 + 2\text{HCl}$ <u>Confirmatory Test Sulphate</u>	white precipitate which is insoluble in dil. HCl.	Presence of sulphate (SO_4^{2-})
1.	To the test soln add lead acetate $\text{M}_2\text{SO}_4 + \text{Pb}(\text{CH}_3\text{COO})_2 \rightarrow \text{PbSO}_4 + 2\text{CH}_3\text{COOH}$	white precipitate	Presence of SO_4^{2-}
2.	To the above precipitate add ammonium acetate solution & heat it. $\text{PbSO}_4 + 2\text{CH}_3\text{COONH}_4 \rightarrow \text{Pb}(\text{CH}_3\text{COO})_2 + (\text{NH}_4)_2\text{SO}_4$	white precipitate dissolves	Presence of SO_4^{2-}

Analysis of cation intergroup separation

- | | | | |
|----|---|-----------------------------|-------------------------|
| 2. | To the above precipitate add ammonium acetate solution & heat it. | white precipitate dissolves | Presence of SO_4^{2-} |
|----|---|-----------------------------|-------------------------|




Analysis of cation intergroup separation

- | | | | |
|----|--|----------------------------------|---|
| 1. | To the test soln add dil. NH_4CO_3 | white precipitate | Absence of 0 group (NH_4^+). |
| 2. | To the test solution add dil. HCl . | No white precipitate. | Absence of I group (Pb^{2+}). |
| 3. | To the test solution add dil. HCl & pass H_2S gas. | No black precipitate. | Absence of II group (Cu^{2+}). |
| 4. | To the test solution add dil. HCl , NH_4Cl & NH_4OH . | No white gelatinous precipitate. | Absence of III group (Al^{3+} , Fe^{2+}). |
| 5. | To the test soln add dil. HCl , NH_4Cl , NH_4OH & pass H_2S gas. | white precipitate | Presence of IV group (Zn^{2+} , Mn^{2+}). |
| 6. | To the test solution add dil. HCl , NH_4Cl , NH_4OH and $(NH_4)_2CO_3$ | No white precipitate. | Absence of V group (Ca^{2+} , Ba^{2+} , Sr^{2+}). |

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
1.	To the test solution add dil. HCl, NH_4Cl , NH_4OH and Na_2HCO_3	no white precipitate	Absence of II group (Hg^{2+})
	Confirmatory test for Zinc ion		
1)	Potassium ferrocyanide solution is added to the test solution	white precipitate	Presence of Zn^{2+}
	$2\text{Zn}^{2+} + \text{K}_4[\text{Fe}(\text{CN})_6] \rightarrow \text{Zn}_2[\text{Fe}(\text{CN})_6] + 4\text{K}^+$ white or bluish white precipitate.		
2)	Ash test The test solution is mixed with few drops of conc. HNO_3 and a little cobalt nitrate solution.	Green tinted ash.	Presence of Zn^{2+} confirmed.

A rolled filter paper dipped in to it is burnt to ash.

the test solution




2) Ash test

The test solution is mixed with few drops of conc MnO_2 and a little cobalt nitrate solution.

Green tinted ash.

Presence of Zn^{2+} confirmed.

A rolled filter paper dipped in to it is burnt to ash.



Result

The given salt contain zinc as cation & sulphate as anion. Therefore the salt is Zinc Sulphate.

~~25/05/26~~

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	<u>Salt 5</u>		
1.	Note the colour of the given salt	Colourless	Absence of Cu^{2+} , Ni^{2+} , Mn^{2+} , Fe^{2+} .
2.	Note the smell of the given salt.	No characteristic smell	Absence of NH_4^+ , CH_3COO^-
3.	Solubility in dil. HCl	Insoluble	Possence of Pb^{2+}
4.	Flame test A pinch of salt is made into paste with conc. HCl & the paste is introduced into the flame by using a glass rod.	No characteristic colour on the flame.	Absence of Ca^{2+} , Ba^{2+} , Sr^{2+} .

4. Flame test

A pinch of salt is made into paste with conc. HCl & the paste is introduced into the flame by using a glass rod.

No characteristic colour on the flame.

Absence of Ca^{2+} , Ba^{2+} , Sr^{2+} .

5. Ash test

The test solution is mixed with few drops of conc. HNO_3 and a little cobalt nitrate solution. A rolled filter paper dipped in it is burnt to ash.

No characteristic colour on the ash.

Absence of Al^{3+} , Zn^{2+} , Mg^{2+} .

Analysis of anion

1) To the test soln add dil. HCl.

No brisk effervescence.

Absence of carbonate.

2) The test solution is fubbed with dil. H_2SO_4 .

No smell of vinegar.

Absence of acetate.

3) To the little of the salt soln about 1ml of conc. H_2SO_4 is added & heated.

No dense white gas.

Absence of chloride.

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
4)	To the test soln add BaCl_2 soln	No white precipitate	Absence of sulphate
5)	A little of the salt solution is heated with conc. H_2SO_4 & a paper ball is added.	Reddish brown fumes.	Presence of nitrate (NO_3^-).
$\text{MNO}_3 + \text{H}_2\text{SO}_4 \longrightarrow \text{MH}_2\text{SO}_4 + \text{HNO}_3$ $4\text{HNO}_3 + \text{C} \longrightarrow 4\text{NO}_2\uparrow + \text{CO}_2 + 2\text{H}_2\text{O}$ (paper ball)			

Confirmatory test

Nitrate

- The test solution is mixed with freshly prepared ferrous sulphate solution. Conc. H_2SO_4 is then added.

A dark brown ring is formed at the junction of two liquids.

Presence of nitrate.

slowly through the sides of the test tube.

(paper ball)

Confirmatory test




Nitrate

1. The test solution is mixed with freshly prepared ferrous sulphate solution. Conc. H_2SO_4 is then added

A dark brown ring is formed at the junction of two liquids.

Presence of nitrate.

slowly through the sides of the test tube.





(Brown ring)

2. Small quantity of solid salt is heated with conc. H_2SO_4 & few copper chips.

Reddish brown fumes

Presence of nitrate confirmed.



(Reddish brown fumes)

Analysis of cation Inter group separation

1. To the test solution add white precipitate dil. Na_2CO_3 .

Absence of O group (NH_4^+)

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
2.	To the test soln add dil. HCl.	White precipitate	Presence of I group (Pb^{2+})
3.	To the test soln add dil. HCl & pass H_2S gas.	No Black precipitate	Absence of II group (Cu^{2+})
4.	To the test solution add dil. HCl, NH_4Cl and NH_4OH & pass H_2S .	No white precipitate	Absence of III group (Al^{3+} , Fe^{2+})
5.	To the test solution add dil. HCl, NH_4Cl & NH_4OH & $(NH_4)_2CO_3$.	No white precipitate.	Absence of IV group (Zn^{2+} , Mn^{2+})
6.	To the test solution add dil. HCl, NH_4Cl & NH_4OH & $(NH_4)_2CO_3$.	No white precipitate	Absence of V group (Ca^{2+} , Ba^{2+} , Sr^{2+})


7.	To the test solution add dil. HCl, NH_4Cl , NH_4OH &	No white precipitate.	Absence of VI group (Mg^{2+})
----	--	-----------------------	-----------------------------------

6. To the test solution add dil. HCl, NH_4Cl & NH_4OH & $(\text{NH}_4)_2\text{CO}_3$	No white precipitate	Absence of V group (Ca^{2+} , Ba^{2+} , Sr^{2+})
--	----------------------	---

7. To the test solution add dil. HCl, NH_4Cl , NH_4OH & Na_2HPO_4	No white precipitate	Absence of VI group (Mg^{2+})
---	----------------------	--

Confirmatory test for lead ion (Pb^{2+})

- | | | |
|---|--------------------|------------------------------|
| 1) The test solution is acidified with acetic & potassium chromate is added | Yellow precipitate | Presence of Pb^{2+} |
|---|--------------------|------------------------------|



- | | | |
|--|--|---|
| 2) Potassium iodide soln is added to the test soln. The above precipitate is boiled with water & cooled. | Yellow precipitate
$\text{Pb}^{2+} + 2\text{KI} \longrightarrow \text{PbI}_2 \downarrow + 2\text{K}^+$
Precipitate dissolves in boiling golden yellow spangles is formed on cooling. | Presence of Pb^{2+}
Presence of Pb^{2+} confirmed. |
|--|--|---|

Result

The given salt contains lead as cation & nitrate as anion. The salt is lead nitrate.

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	Salt 6 Preliminary test		
1.	Note the colour of the given salt.	Colourless.	Absence of Cu^{2+} , Ni^{2+} , Mn^{2+} , Fe^{2+} .
2.	Note the smell of the given salt.	No characteristic smell.	Absence of NH_4^+ and CH_3COO^- .
3.	Solubility in dil. HCl	Soluble	Absence of Pb^{2+} .
4.	Flame test A pinch of salt is made into paste with conc. HCl & the paste is introduced to the flame by using a glass rod.	No characteristic colour on the flame.	Absence of Ca^{2+} , Sr^{2+} , Ba^{2+} .

5) Ash test

Blue tinted

Presence of K^+

introduced to the flame by using a glass rod.

flame.

1) Test

The soln is mixed with few drops of conc HNO_3 & cobalt nitrate. A rolled filter paper dipped & burnt to ash.

Blue tinted ash.

Presence of Al^{3+}

Analysis of anion

1) To the test soln add dil. HCl

No brisk effervescence

Absence of CO_3^{2-}

2) To the test soln is rubbed with dil. H_2SO_4 .

No smell of vinegar

Absence of CH_3COO^-

3) To a little of salt soln about 1ml of conc H_2SO_4 is added & heated.

No dense white gas

Absence of Cl^-

4) A little of the salt soln is heated with conc H_2SO_4 and paper ball is added.

No reddish brown fumes.


Absence of NO_3^-

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
5)	To the test soln add BaCl_2 soln. $\text{M}_2\text{SO}_4 + \text{BaCl}_2 \longrightarrow \text{BaSO}_4 + 2\text{HCl}$	white precipitate which is insoluble in dil HCl.	Presence of sulphate (SO_4^{2-})
	<u>Confirmatory Tests</u> <u>Sulphate</u>		
1)	To the test soln add lead acetate. $\text{M}_2\text{SO}_4 + \text{Pb}(\text{CH}_3\text{COO})_2 \longrightarrow \text{PbSO}_4 + 2\text{CH}_3\text{COOH}$	white precipitate	Presence of SO_4^{2-}
2)	To the above precipitate add ammonium acetate soln and heat it $\text{PbSO}_4 + 2\text{CH}_3\text{COONH}_4 \longrightarrow \text{Pb}(\text{CH}_3\text{COO})_2 + (\text{NH}_4)_2\text{SO}_4$	white precipitate dissolves	Presence of sulphate confirmed.

Analysis of cation Inter group separation

1) To the test soln add white precipitate. Presence of D group


- 2) To the above precipitate add ammonium acetate soln and heat it. white precipitate dissolves. Presence of sulphate confirmed.



Analysis of cation Inter group separation


1) To the test soln add dil Na_2CO_3	white precipitate	Absence of 0 group (NH_4^+)
2) To the test soln add dil. HCl	No white precipitate.	Absence of I group (Pb^{2+})
3) To the test soln add dil. HCl, and pass H_2S gas.	No Black precipitate.	Absence of II group (Cu^{2+})
4) To the test soln add dil HCl, NH_4Cl & NH_4OH .	white gelatinous precipitate.	Presence of III group (Zn^{2+} , Mn^{2+})
5) To the test soln add dil HCl, NH_4Cl , NH_4OH Pass H_2S .	No white precipitate	Absence of IV group (Zn^{2+} , Mn^{2+})
6) To the test soln add dil. HCl, NH_4Cl , NH_4OH & $(\text{NH}_4)_2\text{CO}_3$	No white precipitate	Absence of V group (Ca^{2+} , Ba^{2+} , Sr^{2+})
7) To the test soln add dil HCl, NH_4Cl , NH_4OH & Na_2HCO_3 .	No white precipitate	Absence of VI group (Mg^{2+})

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	<u>Confirmatory test for Aluminium ion (Al^{3+})</u>		
1)	Sodium hydroxide soln. is added in drops to excess to the test soln.	A white precipitate is formed which dissolves in excess of NaOH.	Presence of Al^{3+}
	$Al^{3+} + 3NaOH \longrightarrow Al(OH)_3 + 3Na^+$ $Al(OH)_3 + NaOH \longrightarrow NaAlO_2 + 2H_2O$		
2)	Ash test. The test soln is mixed with a few drops of conc. HNO_3 and a little cobalt nitrate soln. A rolled filter paper is dipped in the soln and burnt to ash.	Blue tinted ash	Presence of Al^{3+} confirmed



Result


of NaOH.



2) Ash test. The test soln is mixed with a few drops of conc HNO_3 and a little cobalt nitrate soln. A rolled filter paper is dipped in the soln and burnt to ash.

Blue tinted ash

Presence of Al^{3+} confirmed



Result

The given salt contain aluminium as cation & sulphate as anion. Therefore the given salt is $\text{Al}_2(\text{SO}_4)_3$.

~~25~~ 25

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	Salt 7 Preliminary Test		
1)	Note the colour of the given salt	Colourless ✓	Absence of Cu^{2+} , Ni^{2+} , Mn^{2+} , Fe^{2+}
2)	Note the smell of the given salt	No characteristic smell	Absence of NH_4^+ , CH_3COO^-
3)	Solubility in dil. HCl	Soluble	Absence of Pb^{2+}
4)	Flame test A pinch of salt is made conc. HCl and the paste is introduced into the flame by using a glass rod.	Pale green flame	Absence of Ba^{2+}
5)	Ash test	No characteristic	Absence of

3) Solubility in dil. HCl	Soluble	Absence of Pb^{2+}
4) Flame test A pinch of salt is made conc. HCl and the paste is introduced into the flame by using a glass rod.	Pale green flame	Presence of Ba^{2+}

5) Ash test The test solution is mixed with few drops of conc. HNO_3 & a little cobalt nitrate soln. A rolled filter paper dipped into it is burnt to ash.	No characteristic colour noted in the ash.	Absence of Al^{3+} , Zn^{2+} , Mg^{2+} .
<u>Analysis of Anions</u>		
1) To the test soln add dil. HCl.	No brisk effervescence	Absence of carbonate (CO_3^{2-}).
2) The test soln is rubbed with dil. H_2SO_4 .	No smell of vinegar	Absence of nitrate (NO_3^-).
3) To a little of the salt soln about 1ml of conc. H_2SO_4 is added & heated.	Dense white gas which fumes with NH_4OH soln taken on a glass rod.	Presence of chloride (Cl^-)

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
	$MCl + H_2SO_4 \longrightarrow MH_2SO_4 + HCl$ $HCl + NH_4OH \longrightarrow NH_4Cl$		
4)	A little of the salt soln is heated with conc H_2SO_4 and paper ball is added	No reddish brown fumes.	Absence of nitrate (NO_3^-)
5)	To the test soln add $BaCl_2$ soln.	No white precipitate	Absence of sulphate (SO_4^{2-})
	<u>Confirmatory Test.</u> <u>Chloride</u>		
1)	The test soln is heated with conc H_2SO_4 and a pinch of a manganese dioxide.	Greenish yellow gas with irritating smell	Presence of chloride
	$2MCl + MnO_2 + 3H_2SO_4 \longrightarrow 2MH_2SO_4 + MnSO_4 + Cl_2 + 2H_2O$		

2) The test soln is mixed white muddy Presence of chloride confirmed

1) The test soln is heated with conc H_2SO_4 and a pinch of a manganese dioxide.	Greenish yellow gas with irritating smell	Presence of chloride
$2MCl + MnO_2 + 3H_2SO_4 \longrightarrow 2HHSO_4 + MnSO_4 + Cl_2 + 2H_2O$		

2) The test soln is mixed with dil HNO_3 and add $AgNO_3$ soln.	white curdy precipitate soluble in excess of NH_4OH .	Presence of chloride confirmed.
$MCl + AgNO_3 \longrightarrow AgCl + MNO_3$ $AgCl + 2NH_4OH \longrightarrow [Ag(NH_3)_2]Cl + 2H_2O$		

Analysis of cations Intergroup Separation

1) To the test soln add dil NH_3 .	white precipitate.	Absence of 0 group (NH_4^+)
2) To the test soln add dil. HCl	No white precipitate	Absence of I group (Pb^{2+})
3) To the test soln add dil. HCl , pass H_2S gas.	No Black precipitate.	Absence of II group (Cu^{2+})
4) To the test soln add dil. HCl , NH_4Cl and NH_4OH	No white gelatinous precipitate	Absence of III group (Al^{3+} , Fe^{3+})

SL. NO.	EXPERIMENT	OBSERVATION	INFERENCE
5)	To the test soln add dil. HCl, NH_4Cl , NH_4OH & pass H_2S .	No white or flesh colored precipitate.	Absence of IV group (Fe^{2+} , Mn^{2+})
6)	To the test soln add dil. HCl, NH_4Cl , NH_4OH & $(\text{NH}_4)_2\text{CO}_3$.	White precipitate.	Absence of V group (Ca^{2+} , Ba^{2+} , Sr^{2+})
7)	To the test soln add dil. HCl, NH_4Cl , NH_4OH & Na_2HPO_4 .	No white precipitate.	Absence of VI group (Mg^{2+})
<u>Confirmatory test for fifth group</u>			
1)	To the test solution acetic acid and potassium chromate soln is added.	Yellow precipitate.	Presence of Ba^{2+} .
	$\text{Ba}^{2+} + \text{K}_2\text{CrO}_4 \longrightarrow$	$\text{BaCrO}_4 + 2\text{K}^+$	

- 2) Ammonium oxalate soln & NH_4OH soln is added to the test soln. The above white precipitate readily dissolved in acetic acid. Presence of Ba^{2+} confirmed.

1) To the test solution acetic acid and potassium chromate soln is added.	Yellow precipitate.	Presence of Ba^{2+} .
$Ba^{2+} + K_2CrO_4 \longrightarrow BaCrO_4 + 2K^+$		

2) Ammonium oxalate soln & NH_4OH soln is added to the test soln. The above precipitate acetic acid is added & heated.	white precipitate readily dissolved in acetic acid. Precipitate readily dissolved.	Presence of Ba^{2+} confirmed. Presence of Ba^{2+} .
$Ba^{2+} + (NH_4)_2C_2O_4 \longrightarrow BaC_2O_4 + 2NH_4^+$		

3) Flame test A pinch of salt is made in the paste with conc. HCl and the paste is introduced in the flame by using a glass rod.	Pale green flame.	Presence of Ba^{2+} .
---	-------------------	-------------------------

Result

The given contains Barium as cation and chloride as anion. Therefore the given salt is $BaCl_2$.

25/05/20

